



NEW

ORBITAL SHAKER



3500.EU.01

- LED DISPLAY
- PLASTIC
- LOW VOLTAGE POWER
- DC MOTOR
- AUDIO ALARM

RPM
50-300 RPM

Max. Capacity
4 kgs

Body Material
Plastic

Glassco Laboratory Equipments Pvt. Ltd.

CE
certified



info@glasscolabs.com | www.glasscolabs.com
+91 171 2891633, 70820-15113, 70820-15114

50
YEARS OF
Celebration
SINCE 1969


ORBITAL SHAKER

Product Features

- Compact design
- Available with Flat platform (Cat. No. 3500.303.01)
- Built in DC brushless motor
- Microprocessor controlled
- User-Friendly
- Application in molecular biology, immunology, chemistry lab
- Low power consumption
- Maximum stirring speed upto 300RPM
- Digital set up of RPM and time with independent LED display
- Audible alarm at time cycle end for safety and ease of unsupervised operation
- Silent and stable operation
- Compatible for use with various capacity conical flasks and different sample containers.

General Optional Accessories

Platform, Flat, Non-Slip Mat



3500.303.01

Platform with clamps, 6 Places (250ml)



3500.303.04

Platform, Double



3500.303.02

Platform with clamps, 12 Places (100ml)



3500.303.05

Platform, Universal, with 3 secured bars



3500.303.03

Platform universal for tubes



3500.303.06

Technical Specification :

Shaker Body Material	: Plastic
Speed Range	: 50-300 RPM
Orbit	: 20 MM
Dimensions of Body	: 284 x 264 x 108 mm
Max. Loading Capacity	: 2.0kgs
Permissible ambient Temp.	: 5 - 45°C
Voltage	: AC 100 - 230V (50/60 Hz)
Timing Range	: 1min - 99Hrs 59Min

Model No.	Orbital Shaker comes	V	Hz
3500.EU.01	with platform, Flat, Non-Slip Mat	230 V	50

Glassco Laboratory Equipments Pvt. Ltd.